# AEROSURE MEDIC<sup>+</sup>

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### Références

1. Padula Ca et al. Inspiratory Muscle Training: Integrative Review. Research and theory for nursing practice 2006; 20:291-304
2. Gosselink et al. Impact of inspiratory muscle training in patients with COPD: what is the evidence? Eur Respir J 2011; 37: 416-425
3. Silva I S et al. Inspiratory muscle training for asthma (Review). 2013 The Cochrane Collaboration. Published by JohnWiley & Sons, Ltd.
4. Bradley JM, Moran F, Greenstone M. Physical training for bronchiectasis (Review). 2002 The Cochrane Collaboration. Published by JohnWiley & Sons, Ltd.
5. Weiner P, Berar-Yanay N, Davidovich A, Magadle R, Weiner M. Specific inspiratory muscle training in patients with mild asthma and high consumption of inhaled beta(2)-agonists. Chest 2000; 117:722-727
6. Lisboa c, Munoz V, Beroiza T, Leiva A, Cruz E. Inspiratory muscle training in chronic airflow limitation: comparison of two different training loads with threshold device. Eur Respir J 1994; 7: 1266-1274
7. Gething AD, Williams M, Davies B. Inspiratory resistive loading improves cycling capacity: a placebo controlled trial. Br J Sports Med 2004; 38: 730-736.
8. Griffiths LA, McConnell AK. The influence of inspiratory and expiratory muscle training upon rowing performance. Eur J Appl Physiol 2007; 99: 457-466.
9. Goosey-Tolfrey V, Foden E, Perret C, Degens H. Effects of inspiratory muscle training on respiratory function and repetitive sprint performance in wheelchair basketball players. Br J Sports Med 2010; 44: 665-668.
10. Magadle R, McConnell AK, Beckerman M, Weiner P. Inspiratory muscle training in pulmonary rehabilitation program in COPD patients. Respir Med 2007; 101: 1500-1505.
11. Konstan M W, Stern R C, Doershuk C F. Efficacy of the flutter device for airway mucus clearance in patients with cystic fibrosis. J Pediatrics 1994; 124: 689-693
12. Endo H, et al. Effects of chest wall vibration on breathlessness during hypercapnic ventilatory response. J Appl Physiol 1998; 84: 1487-1491
13. Fujie T, et al. Effect of chest wall vibration on dyspnea during exercise in chronic obstructive pulmonary disease. Respir Physiol Neurobiol 2002; 130: 305-316
14. App E M, et al. Sputum Rheology Changes in Cystic Fibrosis Lung Disease Following Two Different Types of Physiotherapy: Flutter vs Autogenic Drainage. Chest 1998; 114: 171-177
15. Bhowmik A et al. Improving mucociliary clearance in chronic obstructive pulmonary disease. Respir Med 2009; 103: 496-502
16. Murray MP, Pentland J L, Hill A T. A randomised crossover trial of chest physiotherapy in non cystic fibrosis bronchiectasis. Eur Respir J 2009; 34: 1086-1092
17. Bott J et al. Guidelines for the physiotherapy management of the adult, medical, spontaneously breathing patient. Thorax 2009; 64:i1-i51. Doi:10.1136/thx.2008.110726
18. National Institute of Health and Care Excellence. Chronic obstructive pulmonary disease: Management of chronic obstructive pulmonary disease in adults in primary and secondary care. 2010
19. Sumners D P et al. Increases in inspiratory neural drive in response to rapid oscillating airflow braking forces (vibration). Resp Physiol & Neurobiol 2008; 160: 350-352
20. Sumners D P et al. Respiratory modulation in response to high-frequency airway occlusion delivered during inspiration or expiration. Eur J Appl Physiol 2009. DOI 10.1007/s00421-009-1060-3
21. Morris T, Sumners D P, Green D A. Inspiratory high frequency airway oscillation attenuates resistive loaded dyspnea and modulates respiratory function in young healthy individuals. PLOS ONE 2014. Vol 9, Issue 3

**actegy**<sup>®</sup>  
health powered by science  
inspired by you<sup>™</sup>

### Actegy Ltd.

Reflex, Cain Road, Bracknell, RG12 1HL, Royaume-Uni

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